



Single-use instruments. No risk, all the convenience.

Our Instruments are:

- Quick and easy to set up and use
- Prewrapped in a range of convenient kits and instruments
- Sterile at the point of use
- Convenient to dispose of in sharps container
- Reducing the risk of sharps injuries
- Sharp where needed, ensuring the highest quality of treatment for every procedure
- Implemented to significantly reduce the risk of cross-contamination
- Can be tracked and audited via manufacture and expiry dates on the packaging
- GS1 compliant with removable tracking labels
- TGA, CE, FDA registered
- 5 year expiry date

MDDI Can Increase Revenue by:

- Reducing the number of decontamination cycles per day - minimising costs on consumables, test kits and utilities
- Putting the RDN back chairside to assist the dentist and patient
- Increasing the number of patients seen in each clinical session
- Improving clinical efficiencies allowing for longer treatment plans and procedures
- Improving the patient journey

STERILE, SINGLE-USE DIAGNOSTIC INSTRUMENTS

MDDI kits contain high quality, precision made instruments. Our range of probes and our college tweezers are manufactured from stainless steel and all our handles are moulded with recyclable ABS plastic.



Front Face Mirror with Percussion Tip

SI053ST-1

180 mirrors per dispenser box



Tweezers

SI081ST-1

200 tweezers per dispenser box



London College Tweezers

SI182ST-1

200 tweezers per dispenser box



Periodontal 2 Piece Examination Kit

SI552ST-1

100 kits per dispenser box

1 x Front Face Mirror with Percussion Tip
1 x Perio Probe/Sickle



Periodontal 3 Piece Examination Kit

SI551ST-1

100 kits per dispenser box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x Tweezers



Periodontal 3 Piece Examination Kit with CPITN-C Probe

SI951ST-1

100 kits per dispenser box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No.8 Standard Probe/CPITN-C Probe
3.5/5.5/8.5/11.5
- 1 x Tweezers



Premium Periodontal 3 Piece Examination Kit

SI631ST-1

100 kits per dispenser box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x London College Tweezers



Periodontal 2 Piece Examination Kit with CPITN-C Probe

SI952ST-1

100 kits per dispenser box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No.8 Standard Probe/CPITN-C Probe
3.5/5.5/8.5/11.5



Complete Periodontal Examination Kit

SI550ST-1

25 kits per inner box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x Tweezers
- 2 x Face Mask with Ear Loops
- 1 x Patient Bib with Ties
- 2 x Cotton Roll No.2
- 1 x Tissue Wipe
- 1 x Tray (19.5cm x 12.5cm x 2cm)



Complete Periodontal Examination Kit with CPITN-C Probe

SI950ST-1

25 kits per inner box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No.8 Standard Probe/CPITN-C Probe 3.5/5.5/8.5/11.5
- 1 x Tweezers
- 2 x Face Mask with Ear Loops
- 1 x Patient Bib with Ties
- 2 x Cotton Roll No.2
- 1 x Tissue Wipe
- 1 x Tray (19.5cm x 12.5cm x 2cm)



Suture Removal Kit

SI770ST-1

25 kits per inner box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No.8 Standard Probe/CPITN-C Probe 3.5/5.5/8.5/11.5
- 2 x Face Mask with Ear Loops
- 2 x Non-Woven Swab 4ply (5cm x 5cm)
- 1 x Iris Scissors Curved 11.5cm
- 1 x Iris Toothed Forceps 10cm
- 1 x Tray (19.5cm x 12.5cm x 2cm)
- 1 x Tweezers

STERILE, SINGLE-USE RESTORATIVE KITS

2 Piece Restorative Mixing Kit

SI622ST-1

150 kits per dispenser box

1 x Spatula DE
1 x Dycal Applicator



3 Piece Composite Restorative Kit

SI630ST-1

100 kits per dispenser box

1 x Flat Plastic DE
1 x Ball Burnisher DE
1 x Medium Spoon Excavator DE (1.8mm)



5 Piece Amalgam Restorative Kit

SI650ST-1

60 kits per dispenser box

1 x Flat Plastic DE
1 x Ball Burnisher DE
1 x Medium Spoon Excavator DE (1.8mm)
1 x Carver DE
1 x Amalgam Plugger Large/Small



"From a time, organisation and environmental cleaning perspective, the disposable instruments are fantastic. Quicker changeovers mean more patients seen and less staff time utilised in the sterilisation bay, which is less handling of sharp instruments too."

- Dental Assistant - Oral Health Outreach

6 Piece Composite Restorative Kit

SI660ST-1

50 kits per dispenser box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Small Spoon Excavator DE (1.2mm)
- 1 x Medium Spoon Excavator DE (1.8mm)
- 1 x PFI49 Burnisher
- 1 x 1/2 Hollenback DE



6 Piece Comprehensive Restorative Kit

SI661ST-1

50 kits per dispenser box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Medium Spoon Excavator DE (1.8mm)
- 1 x Carver DE
- 1 x Ladmore Plugger DE
- 1 x Amalgam Plugger Large/Small



MDDI's restorative and hygienist instruments have an ergonomic, ABS handle moulded around a stainless steel core, ensuring the instruments display strength and rigidity when required.

