

CLINICAL CONSUMABLES CATALOGUE

www.mddiglobal.com



CONTENTS

SURGICAL DRAPES	3
.....	
TUBING	4
.....	
FACE MASKS AND OTHER CONSUMABLES	5
.....	
IMPRESSION TRAYS	6
.....	

MDDI is constantly expanding our product range so if you are having difficulty sourcing a sterile single-use product to suit your needs, please do not hesitate to contact us.

sales@mddiglobal.com

In response to critical demand MDDI has developed a cost effective and innovative range of single-use clinical consumables suitable for use within both dental and medical environments.

SURGICAL DRAPES



Surgical Gown Bundle

PA073ST-1

30 gown bundles per carton

1 x Surgeon Gown (125cm x 155cm)

60gsm AAMI Level 2

2 x Hand Towel (40cm x 40cm)



Sterile Surgical Drape Impervious (90cm x 120cm)

SD011ST-1

18 drapes per box

8 boxes per carton



Sterile Surgical Drape Impervious (75cm x 90cm)

SD013ST-1

25 drapes per box

8 boxes per carton




Sterile 3/4 Patient Drape/U Drape with Adhesive + Tube Holders (100cm x 150cm)

SD016ST-1

10 drapes per box

8 boxes per carton





**Aspirator Kit
CC072ST-1**

8 kits per box | 8 boxes per carton


- 1 x Suction Tube with Connectors (3m)
- 1 x Fine Aspiration Handle and Tip with Control Vent (2.8mmID)

Irrigation Tubing - Set A

CC101ST-1

20 sets per box | 8 boxes per carton

Compatible with NSK – Surgic XT/Surgic XT PLUS/Variosurg, Biomet 3i – Osseocision, KYOCERA – Implantor Neo. ASEPTICO – AEU 7000/AEU 6000 motors.



Irrigation Tubing - Set B

CC102ST-1

20 sets per box | 8 boxes per carton


Compatible with Nobel Biocare80, W&H Implantmed, Implant Med new design 2005, Implantmed SI 923, Astratech Implant SI 923, Straumann, Osseoset 100, Osseoset 200, Osseocare.



**Surgical Face Mask
with Elasticated
Ear Loops - Level 3**

PA014NS-1

50 masks per box
40 boxes per carton



Our face masks are CE certified conforming to AS 4381:2015

- Non-sterile • Latex free • Earloops • 3-ply non-woven fabric with 25gsm filtration
- Fluid resistance 160mmHg • Delta-P <5.0 • PFE >99% @ 0.1µm • BFE >98%

BARRIERS



**Drill Sleeve Elasticised Clear on Card
(120cm x 7.5cm) Pkt1**

CC126ST-1

50 individual sleeves packaged as 50 packs of 1
8 boxes per carton



**Drill Sleeve Elasticised Clear on Card
(120cm x 7.5cm) Pkt2**

CC026ST-1

60 individual sleeves packaged as 30 packs of 2
8 boxes per carton



**Quick Drill Sleeve with Applicator
(150cm x 9.5cm) Pkt1**

CC175ST-1

24 individual sleeves packaged as 24 packs of 1
8 boxes per carton



**Quick Drill Sleeve with Applicator Ring
(150cm x 9.5cm) Pkt2**

CC075ST-1

24 individual sleeves packaged as 12 packs of 2
8 boxes per carton



Adhesive Film Blue (20cm x 20cm) Pkt 3

CC022ST-1

50 packs per box | 6 boxes per carton



Single-use plastic impression trays with integrated handle. Constructed from rigid, recyclable ABS for accurate, comfortable and distortion free impressions.

Non-Sterile. 25 trays per box | 4 boxes per carton



Impression Tray Small Upper
IT02SMU-1



Impression Tray Medium Upper
IT04MDU-1



Impression Tray Large Upper
IT06LGU-1



Impression Tray Extra Large Upper
IT08XLU-1



Impression Tray Small Lower
IT01SML-1



Impression Tray Medium Lower
IT03MDL-1



Impression Tray Large Lower
IT05LGL-1



Impression Tray Extra Large Lower
IT07XLL-1



These products can be purchased via our website
or directly through our sales team.

View our online store at www.mddiglobal.com



Our valued distributors

We are constantly looking for new ways to collaborate,
please do not hesitate to get in touch.

MDDI - HEAD OFFICE

E: sales@mddiglobal.com | P: +61 (0)7 5453 4433

Unit 1, 24 Hoopers Road, Kunda Park Queensland 4556 AUSTRALIA